


we desire every patients with diabetes can be healthy



Insulin pen needle contrast before and after using



Before using After once use After twice use After sixth use

Type Specification

| Type | Specification |
|----------------|---------------|
| 33G NEW | 0.20x4mm |
| 32G NEW | 0.23x4mm |
| 31G | 0.25x5mm |
| 31G | 0.25x8mm |
| 30G | 0.30x5mm |
| 30G | 0.30x8mm |
| 29G | 0.33x12mm |

Remark: normal size in above table, we can also make special size as client' s request



KINDLY ENTERPRISE

Zhejiang Kindly Medical Devices Co., Ltd

Add:No.758,5th Binhai Road,Binhai Industrial Park,Longwan


District Wenzhou,Zhejiang,China.

Postcode: 325025

Tel: 86-577-86862296 86862297

Fax: 86-577-86856288

www.kdlchina.com



Disposable
Insulin Pen Needle



Insulin Pen Needle Features

Commonality Suit any kind of insulin pen and easy for medical workers and patients to use.

All specifications The product is complete in specification and can satisfy different injection requirement of patients.

Comfortableness Reduce painful feeling of puncture and make patients feel comfortable during injection.

Professional KDL is one of the largest medical needle manufacturers in China.

Product Features



Safety design and easy for using

The strength design in needle set make it more convenient and safe to use

High compatibility and cover almost all branches

The product has high compatibility. The standard design of internal threaded ensure it can stably match all branches insulin pen.

Best penetration make the injection more comfortable

KDL insulin pen needle has a series of advanced and mature needle tip technique and silication technique. The sharp needle tip make puncture more easily and let patients feel more comfortable during injection.

High Matching Degree



NovoPen® 4



NovoPen® 3



NovoPen® Junior



NovoRapid® FlexPen®



BD: Puregon Pen®



Sanofi: Aventis Lantus



Sanofi: OptiPen®



Sanofi: TactiPen®



Owen Mumford:AutoPen®



Owen Mumford:AutoPen® 2



Lily:HumaPen®



Lily:HumaPen®



Wanbang:Wangbang Pen 2




Tong hua dong bao:Dong Bao Pen®



Tong hua dong bao:Gansulin Pen®


The screw thread design of KDL insulin pen needle meets the standards regulated by ISO11608-2, which is compatible with other standard ones domestically produced or imported.

More Comfortable, More Safety



Needle Structure Parsing

Sterile injection guarantee with double protector



33G Superfine Short Insulin Pen Needle

4mm needle is 50% shorter than traditional 8mm needle



Reduce the risk of injecting to muscle layer directly and make the injection more safely.

Don't need to pinch the skin part during injection which make it more convenient and easy in operation.

33G insulin pen needle will create more comfortable injection feeling to patient.

33G, a shorter and finer insulin pen needle!!

